

Ultraclean Sonozyme 3 Technical Data Sheet



Description

Ultrawave's Sonozyme 3 is a pH neutral tri-enzymatic detergent for cleaning and pre-disinfection of surgical instruments and medical devices. As well as containing enzymes for breaking down organic contaminant.

Applications

For use within (but not limited to) healthcare, dental, podiatry, tattooist and veterinary sectors.

Directions for Use

Add cold or tepid water (not exceeding 45°C) to the bath. Add Sonozyme 3 to obtain a dilution of 5mL per Litre, that's 40mL in 8Litres of water. Operate the Ultrasonic bath as directed in the user guide. After processing, rinse the cleaned devices with water. The fluid in the bath should be changed and re-dosed at the start of each clinical session (twice a day) or more frequently if visibly contaminated. Please refer to the product safety data sheet before use.

Additional note: The solution may turn cloudy when heated, this is normal and does not effect the cleaning performance.

Technical Data

These figures are typical and should not be interpreted as specifications, refer to the product safety data sheet for details.

Characteristic	Value/Description
Appearance	Blue liquid
Odour	Characteristic
pH (@ 1%)	7.8 (neutral)
Relative Density at 20°C	1.031

Supplied as 3 x 2Litre Bottles (SKU: F50026).

Disposal

Process effluent can normally be discharged to foul sewer (subject to consent limits). Empty cleaned containers can be recycled where facilities exist.

Storage

Store in a cool, dry place protected from frost and away from acids and strong oxidising agents. Store upright in original containers. Recommended storage temperature 5-25°C. Shelf-life up to 2years from date of manufacture.

Specifications are subject to change for continual improvements