

Ultraclean Disinfect+ Technical Data Sheet



Description

Ultrawave's Disinfect+ is a pH neutral 2 in 1 solution that cleans and disinfects in one. Most disinfectant fluids require you to pre-clean tools, instruments and surfaces before use. Disinfect+ gives you all the power of a disinfectant with added ease of use. As you don't have to pre-clean, you can just use the Disinfect+ in your Ultrasonic bath and your items will be clean and disinfected in as little as 6 minutes*. Independently tested as Bactericidal, Mycobactericidal/Tuberculocidal, Fungicidal/Yeasticidal and effective against envelope viruses (including Coronaviruses). *Requires 15 minutes for Mycobactericidal/Tuberculocidal

Applications

For use in all markets/sectors that require disinfection.

Directions for Use

Dilute with water to give a 5% (50mL per Litre) solution before use with water. Ensure thorough mixing of the solutions. Follow the device manufacturer's instructions for cleaning. Remove gross soiling before disinfection. Apply by immersion or using ultrasonic agitation. Operate the Ultrasonic bath as directed in the user guide. After processing, rinse the cleaned devices with water. It is recommended that the fluid in the bath should be changed and re-dosed twice a day (or more frequently if visibly contaminated). Do not apply to hot surfaces or live electrical equipment. This product is not intended for use on invasive medical devices. Please refer to the product safety data sheet before use. Use only as directed.

Technical Data

These figures are typical and should not be interpreted as specifications, refer to the product safety data sheet for details

Characteristic	Value/Description
Appearance	Clear, colourless liquid
Odour	Virtually odourless
pH (@ 1%)	7.5 (neutral)
Relative Density at 20°C	0.997

Supplied as 2 x 5Litre Bottles (SKU: F50036).

Disposal

Process effluent can normally be discharged to foul sewer (subject to consent limits). Empty cleaned containers can be recycled where facilities exist.

Storage

Store in a cool, dry place protected from frost and away from acids and strong oxidising agents. Store upright in original containers. Recommended storage temperature 5-25°C. Shelf-life up to 2years from date of manufacture.

Specifications are subject to change for continual improvements