

QS5

4.5Litre Ultrasonic Bath

Ultrawave's **Q-Series** offers outstanding performance for use in medical, scientific, laboratory, chiropody, veterinary and many other professional environments. **Q-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

Incorporating Ultrawave's Frequency Leap Technology, homogenous cavitation activity is generated throughout the ultrasonic bath reducing dead spots and standing waves. This provides a more consistent scrubbing action for excellent results.



The **QS5** is a 4.5Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode, fluid level sensor and a memory card port, the **QS5** offers repeatable and traceable cleaning performance with consistent cleaning results every time.

Specifications are subject to change for continual improvements

QS5

4.5Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Digital LCD
Overall Dimensions	340 x 325 x 265mm (WxDxH)
Tank Dimensions	300 x 150 x 150mm (WxDxH)
Internal Basket Dimensions	265 x 115 x 110mm (WxDxH)
Maximum Capacity	5.0Litres
Working Capacity	4.5Litres
Ultrasonic Power	100W
Ultrasonic Power (per litre)	22.2W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 70°C
Heater Power	150W
Timer	0 - 99 minutes
Mains Connection	230VAC
Drain Outlet	1/2" BSP valve (located at rear of machine)
Lid	ABS Plastic
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	6.5kg
Net Weight	5.5kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)