



Ultrawave


Hygea Valet Instrument Management System

The Hygea Valet Instrument management system securely holds dental instruments in place throughout the decontamination process. Incorporating an open design to aid the ultrasonic cleaning process, the Hygea Valet cassettes significantly reduce the potential for sharp injuries and allow easy transfer throughout the surgery.



Incorporating a strong stainless steel frame and removable silicone strip inserts, the Hygea Valet cassettes are available in 3 sizes for holding 5, 10 and 15 dental examination instruments.

These systems will ensure that instruments are protected and stored throughout the whole decontamination cycle as the cassettes can also be used in washer disinfectors and autoclaves.



Hygea Valet Instrument Management System

- Self locking cassettes securely hold instruments in place
- Minimises risk of sharps injuries
- Open cassettes design promotes maximum ultrasonic activity
- Sizes designed for 5, 10 and 15 instruments
- Tailored design specific to the Hygea 2 basket carrier and allows to clean up to 45 standard dental instruments in one cycle.
- Suitable for use in ultrasonic cleaners, washer disinfectors and autoclaves

Hygea Valet Instrument Management System

5 way

10 way

15 way

Dimensions (L x W x H)

180 x 95 x 27

180 x 136 x 27

180 x 215 x 27

Eastgate Business Park,
Wentloog Avenue,
Cardiff, CF3 2EY, UK

Tel: +44 (0) 845 330 4236
Fax: +44 (0) 845 330 4231
admin@ultrawave.co.uk
www.ultrawave.co.uk

