

# QS18

## 17.5Litre Ultrasonic Bath

Ultrawave's **Q-Series** offers outstanding performance for use in medical, scientific, laboratory, chiropody, veterinary and many other professional environments. **Q-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

Incorporating Ultrawave's Frequency Leap Technology, homogenous cavitation activity is generated throughout the ultrasonic bath reducing dead spots and standing waves. This provides a more consistent scrubbing action for excellent results.



The **QS18** is a 17.5Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode, fluid level sensor and a memory card port, the **QS18** offers repeatable and traceable cleaning performance with consistent cleaning results every time.

Specifications are subject to change for continual improvements

# QS18

## 17.5Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Digital LCD
Overall Dimensions	375 x 490 x 365mm (WxDxH)
Tank Dimensions	330 x 300 x 200mm (WxDxH)
Internal Basket Dimensions	295 x 265 x 160mm (WxDxH)
Maximum Capacity	18.5Litres
Working Capacity	17.5Litres
Ultrasonic Power	300W
Ultrasonic Power (per litre)	17W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 70°C
Heater Power	450W
Timer	0 - 99 minutes
Mains Connection	230VAC
Drain Outlet	1/2" BSP valve (located at rear of machine)
Lid	ABS Plastic
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	16.5kg
Net Weight	15kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)