

# Qi300

## 17Litre Ultrasonic Bath

Ultrawave's **Qi-Series** are robust and reliable, delivering outstanding cleaning performance in a wide range of industries. From precision, pharmaceutical and medical device manufacturing to heavy duty cleaning in automotive and engineering processing. **Qi-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

The 2mm thick stainless steel tank provides longevity and durability whilst excellent ultrasonic activity is generated via Ultrawave's *Frequency Leap* technology.

*Frequency Leap* technology reduces dead spots and promotes homogenous cavitation activity within the tank. This provides a more consistent scrubbing action for highly effective cleaning results.



The **Qi-300** is a 17Litre ultrasonic cleaner with a menu driven, LCD control panel. This allows simple adjustment of cycle time (0 to 99mins or infinity mode), temperature (ambient to 80°C), ultrasonic power and frequency mode for precise and accurate processing.

Incorporating features such as degas mode and a fluid level sensor, the

**Qi-Series** is reliable and durable, giving you consistent cleaning results every time.



Specifications are subject to change for continual improvements

# Qi300

## 17Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	316 Stainless Steel
Control	Digital LCD
Overall Dimensions	350 x 380 x 395mm (WxDxH)
Tank Dimensions	305 x 275 x 200mm (WxDxH)
Internal Basket Dimensions	275 x 250 x 125mm (WxDxH)
Maximum Capacity	17Litres
Working Capacity	14.5Litres
Ultrasonic Power	300W
Ultrasonic Power (per litre)	20W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 80°C
Heater Power	750W
Timer	0 - 99 minutes or Continuous
Mains Connection	230VAC
Drain Outlet	Lockable 1/2" BSP valve (located at rear of machine)
Lid	304 Stainless Steel
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	18.8kg
Net Weight	16.3kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)