

hygea® 6427 Accessory Kit

The hygea 6427 accessory kit (SKU S40031) includes a selection of tubing and connectors, that can be used for cleaning a wide range of both laparoscopic and standard instruments.



SKU	Description	Qty
C40215	Luerlok insert, no valve, 3/16in barb	12
C40216	Luerlok insert, no valve, 1/4in barb	12
C40217	Luerlok insert, no valve, 1/8in barb	12
C40221	Male Hose Coupling Straight, 3mm	6
C40223	Female luer to 4mm hosetail	12
C40224	Male luer to 5mm hosetail	12
C40225	Luer threaded female to 4mm hosetail	3
C40226	Luer threaded male to 5mm hosetail	6
C40232	Male Hose Coupling Straight, 4mm	6
C40233	Male Hose Coupling Straight, 5mm	6
C40227	Silicone Tube, 4.0mm ID, 2mm wall	12 80mm lengths
C40228	Silicone Tube, 5.0mm ID, 1.5mm wall	12 100mm lengths
C40231	Silicone Tube, 6.0mm ID, 2mm wall	6 120mm lengths
C40237	Silicone Tube, 3.0mm ID, 2mm wall	12 100mm lengths

Specifications are subject to change for continual improvements