

Qi200

12Litre Ultrasonic Bath

Ultrawave's **Qi-Series** are robust and reliable, delivering outstanding cleaning performance in a wide range of industries. From precision, pharmaceutical and medical device manufacturing to heavy duty cleaning in automotive and engineering processing. **Qi-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

The 2mm thick stainless steel tank provides longevity and durability whilst excellent ultrasonic activity is generated via Ultrawave's *Frequency Leap* technology.

Frequency Leap technology reduces dead spots and promotes homogenous cavitation activity within the tank. This provides a more consistent scrubbing action for highly effective cleaning results.



The **Qi-200** is a 12Litre ultrasonic cleaner with a menu driven, LCD control panel. This allows simple adjustment of cycle time (0 to 99mins or infinity mode), temperature (ambient to 80°C), ultrasonic power and frequency mode for precise and accurate processing.

Incorporating features such as degas mode and a fluid level sensor, the

Qi-Series is reliable and durable, giving you consistent cleaning results every time.

Specifications are subject to change for continual improvements

Qi200

12Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	316 Stainless Steel
Control	Digital LCD
Overall Dimensions	320 x 320 x 395mm (WxDxH)
Tank Dimensions	275 x 217 x 200mm (WxDxH)
Internal Basket Dimensions	245 x 195 x 125mm (WxDxH)
Maximum Capacity	12Litres
Working Capacity	10.5Litres
Ultrasonic Power	200W
Ultrasonic Power (per litre)	19W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 80°C
Heater Power	600W
Timer	0 - 99 minutes or Continuous
Mains Connection	230VAC
Drain Outlet	Lockable 1/2" BSP valve (located at rear of machine)
Lid	304 Stainless Steel
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	15.4kg
Net Weight	13kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)