

Sanitiser 13

Ultrasonic Bath

Ultrawave's **Sanitiser-Series** offers outstanding performance for use in medical, healthcare, scientific, laboratory, chiropody, veterinary and many other professional environments. **Sanitiser-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

Cleans and disinfects, when used in conjunction with Ultraclean Disinfect+.

Effective against enveloped viruses (including Coronaviruses), bactericidal, mycobactericidal/tuberculocidal and fungicidal/yeastocidal.



The **Sanitiser 13** is a 12.5Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode and fluid level sensor, the **Sanitiser 13** offers repeatable cleaning performance with consistent cleaning results every time.

Specifications are subject to change for continual improvements

Sanitiser 13

Ultrasonic Bath

Construction Case	Powder coated 304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Digital LCD
Overall Dimensions	365 x 415 x 265mm (WxDxH)
Tank Dimensions	330 x 300 x 150mm (WxDxH)
Internal Basket Dimensions	295 x 265 x 110mm (WxDxH)
Maximum Capacity	13Litres
Working Capacity	12.5Litres
Ultrasonic Power	200W
Ultrasonic Power (per litre)	16W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 70°C
Heater Power	300W
Timer	0 - 99 minutes
Mains Connection	230VAC
Drain Outlet	1/2" BSP valve (located at rear of machine)
Lid	ABS Plastic
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	8.8kg
Net Weight	7.8kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)