

# Sanitiser 5

## Ultrasonic Bath

Ultrawave's **Sanitiser-Series** offers outstanding performance for use in medical, healthcare, scientific, laboratory, chiropody, veterinary and many other professional environments. **Sanitiser-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

Cleans and disinfects, when used in conjunction with Ultraclean Disinfect+.

Effective against enveloped viruses (including Coronaviruses), bactericidal, mycobactericidal/tuberculocidal and fungicidal/yeastocidal.



The **Sanitiser 5** is a 4.5Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode and fluid level sensor, the **Sanitiser 5** offers repeatable cleaning performance with consistent cleaning results every time.

Specifications are subject to change for continual improvements

# Sanitiser 5

## Ultrasonic Bath

|  |   |
|--|---|
| Construction Case                              | Powder coated 304 Stainless Steel   |
| Construction Tank                              | 304 Stainless Steel   |
| Control  | Digital LCD   |
| Overall Dimensions                             | 335 x 265 x 265mm (WxDxH)   |
| Tank Dimensions                                | 300 x 150 x 150mm (WxDxH)   |
| Internal Basket Dimensions                     | 265 x 115 x 110mm (WxDxH)   |
| Maximum Capacity                               | 5.0Litres   |
| Working Capacity                               | 4.5Litres   |
| Ultrasonic Power                               | 100W  |
| Ultrasonic Power (per litre)                   | 22.2W/Litre   |
| Operating Frequency                            | 32 - 38kHz <i>Frequency Leap</i>  |
| Maximum Heating Capacity                       | Ambient - 70°C  |
| Heater Power                                   | 150W  |
| Timer  | 0 - 99 minutes  |
| Mains Connection                               | 230VAC  |
| Drain Outlet                                   | 1/2" BSP valve (located at rear of machine)   |
| Lid  | ABS Plastic   |
| Basket   | 304 Stainless Steel   |
| Gross Weight<br>(inc. lid, basket & packaging) | 6.5kg   |
| Net Weight                                     | 5.5kg   |
| Warranty                                       | 1 Year as standard - extend to 3<br>Years if registered online with us<br>(within 30days of purchase) |