

Sanitiser 9

Ultrasonic Bath

Ultrawave's **Sanitiser-Series** offers outstanding performance for use in medical, healthcare, scientific, laboratory, chiropody, veterinary and many other professional environments. **Sanitiser-Series** Ultrasonic Cleaners are Class I Medical Devices, compliant to MDR (EU) 2017/745.

Cleans and disinfects, when used in conjunction with Ultraclean Disinfect+.

Effective against enveloped viruses (including Coronaviruses), bactericidal, mycobactericidal/tuberculocidal and fungicidal/yeastocidal.



The **Sanitiser 9** is a 9.5Litre heated, digital ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 70°C), ultrasonic power and frequency mode for precise and accurate processing.

Combining features such as degas mode and fluid level sensor, the **Sanitiser 9** offers repeatable cleaning performance with consistent cleaning results every time.



Specifications are subject to change for continual improvements

Sanitiser 9

Ultrasonic Bath

Construction Case	Powder coated 304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Digital LCD
Overall Dimensions	340 x 265 x 265mm (WxDxH)
Tank Dimensions	300 x 240 x 150mm (WxDxH)
Internal Basket Dimensions	265 x 210 x 110mm (WxDxH)
Maximum Capacity	10Litres
Working Capacity	9.5Litres
Ultrasonic Power	100W
Ultrasonic Power (per litre)	10.5W/Litre
Operating Frequency	32 - 38kHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 70°C
Heater Power	150W
Timer	0 - 99 minutes
Mains Connection	230VAC
Drain Outlet	1/2" BSP valve (located at rear of machine)
Lid	ABS Plastic
Basket	304 Stainless Steel
Gross Weight (inc. lid, basket & packaging)	7.3kg
Net Weight	6.3kg
Warranty	1 Year as standard - extend to 3 Years if registered online with us (within 30days of purchase)